

## **JAX: mutation - Truseq\_JAX,**

### **Molecular Methods Description:**

The variant data (point mutations and indels) has been analyzed from next-generation sequencing using The Illumina TruSeq Amplicon Cancer Panel covering 48 cancer-related genes.

### **Analysis Description:**

- The analysis of the sequencing output uses the Xenome tool to remove contaminating mouse sequences before alignment and variant calling. \* BWA, GATK, and SnpEff are utilized for alignment (GRCh38 human reference), variant calling and annotation. \* The high quality variants are further filtered to remove likely false positives, variants with low and modifier impact based on SnpEff annotations, and putative germline variants. \* More details on the analysis can be found [here](#)